方法 1	方法2
Method 1	Method 2
1.用含 10% 乙醇、0.5% 乙酸的固定液固定 2 次, 6 min/次	1.在含 10%乙醇、0.5%乙酸的 0.2%硝酸银中染色 5min
Place gels in fixing solution including 10% ethanol and 0.5%	Place gels in solution including 10% ethanol, 0.5% acetic acid
acetic acid, and shake gently for 6 min, twice	and 0.2% silver nitrate , and shake gently for 5 min
2.0.2%AgNO ₃ 溶液渗透 10-12min	2.蒸馏水中漂洗 30s
Replace the solution with 0.2% silver nitrate aqueous solution,	Rinse gels briefly with d H ₂ O for 30 s
and shake gently for 10-12 min.	
3.dH ₂ O 清洗 2 次,30s	3.在含 0.5%甲醛的 3.0%氢氧化钠中显色
Rinse gels briefly with d H ₂ O for 30 s, twice	Apply developing solution, including 3.0% Sodium hydroxide,
	0.5% formaldehyde, and shake gently
4.0.002%硫代硫酸钠溶液置大约 30s, 直到溶液变黄, 倒掉	4.条带清晰后,用蒸馏水冲洗数次,终止显色反应
Apply 0.002% sodium thiosulfate aqueous solution, and shake	After the bands appear with desirable intensity, Wash gel in
gently for 30 s. After the solution turned yellow, remove the	distilled water several times to stop the reaction
solution	
5.用含 NaOH 1.5%, 甲醛 0.4%的显色液显色	
Apply developing solution, including 1.5% Sodium hydroxide,	
0.4% formaldehyde, and shake gently	
6.条带清晰后,用蒸馏水冲洗数次,终止显色反应	
After the bands appear with desirable intensity, Wash gel in	
distilled water several times to stop the reaction	
<u> </u>	