

表 1 SRAP 引物序列

Table 1 Primer sequences used for SRAPA anylysis

编号	正向引物	编号	反向引物
Code	Forward primers	Code	Reverse primers
Me1	5'-TGAGTCCAAACCGGATA-3'	Em1	5'-GACTGCGTACGAATTAAT-3'
Me2	5'-TGAGTCCAAACCGGAGC-3'	Em2	5'-GACTGCGTACGAATTTGC-3'
Me3	5'-TGAGTCCAAACCGGAAT-3'	Em3	5'-GACTGCGTACGAATTGAC-3'
Me4	5'-TGAGTCCAAACCGGACC-3'	Em4	5'-GACTGCGTACGAATTTGA-3'
Me5	5'-TGAGTCCAAACCGGAAG-3'	Em5	5'-GACTGCGTACGAATTAAC-3'
Me6	5'-TGAGTCCAAACCGGTAG-3'	Em6	5'-GACTGCGTACGAATTGCA-3'
Me7	5'-TGAGTCCAAACCGGTTG-3'	Em7	5'-GACTGCGTACGAATTATG-3'
Me8	5'-TGAGTCCAAACCGGTGT-3'	Em8	5'-GACTGCGTACGAATTAGC-3'
		Em9	5'-GACTGCGTACGAATTACG-3'
		Em10	5'-GACTGCGTACGAATTTAG-3'