

表 4 ‘Fiesta’和‘Altimo’叶柄外植体的 3 种消毒方法

Table 4 Three method of disinfection for ‘Fiesta’ and ‘Altimo’

组合	消毒方法
Combination	Method of disinfection
A	70%酒精 10 s; 1.5%次氯酸钠浸泡 10 min; 无菌水洗 4 次
	70% Ethanol 10 s; 1.5% NaClO 10 min; Sterile water 4 times
B	50 mg/L 制霉菌素浸泡 20 min; 70%酒精 30 s; 1.5%次氯酸钠浸泡 10 min; 无菌水洗 4 次; 抗氧化剂(150 mg/L 柠檬酸, 100 mg/L 维生素, 200 mg/L 半胱氨酸)处理 10 min
	50 mg/L Nystatin 20 min; 70% Ethanol 30 s; 1.5% NaClO 10 min; Sterile water 4 times; Antioxidant (150 mg/L Citric acid, 100 mg/L Vitamine, 200 mg/L Cysteine) 10 min
C	50 mg/L 制霉菌素浸泡 20 min; 70%酒精 30 s; 0.1%升汞浸泡 5 min; 无菌水洗 4 次
	50 mg/L Nystatin 20 min; 70% Ethanol 30 s; 0.1% Corrosive sublimate 5 min; Sterile water 4 times