

表5用于检测4种原虫的引物序列

Table 5 The primers used to detect four protozoans

原虫 Protozoans	引物 Primers	引物序列(5'-3') Primers sequences (5'-3')	扩增长度(bp) Length of PCR products (bp)
尼氏单孢子虫 Haplosporidium nelsoni	MSX-A	CGACTTTGGCATTAGGTTTCAGACC	573
派琴虫 Perkinsus	MSX-B	ATGTGTTGGTGACGCTAACCG	
折光马尔太虫 Marteilia refringens	PerkITS-85	CCGCTTTGTTTGGATCCC	703
包拉米虫 Bonamia	PerkITS-750	ACATCAGGCCTTCTAATGATG	
	Marteiliosis 1	CCGCACACGTTCTTCACTCC	411
	Marteiliosis 2	CTCGCGAGTTTCGACAGACG	
	CF	CGGGGGCATAATTCAGGAAC	760
	CR	CCATCTGCTGGAGACACAG	