

表 2 Tail-PCR 反应程序

Table 2 Procedure of Tail-PCR

反应	程序号	循环数	反应条件
Reaction	File No.	Cycle No.	Thermal condition
第一轮	1	1	92°C 3 min; 95°C 1 min
1 <sup>st</sup> amplification	2	5	94°C 30 s; 60°C 1 min; 72°C 2 min 30s
	3	1	94°C 30 s; 28°C 2 min, ramping to 72°C 0.3°C/s; 72°C 2 min
	4	10	94°C 30 s; 42°C 1 min; 72°C 2 min 30 s
	5	1	94°C 30 s; 60°C 1 min; 72°C 2 min 30 s
	6	14	94°C 30 s; 60°C 1 min; 72°C 2 min 30 s
			94°C 30 s; 42°C 1 min; 72°C 2 min 30 s
	7	1	72°C 5 min
第二轮	1	1	92°C 3 min; 95°C 1 min
2 <sup>nd</sup> amplification	8	15	94°C 30 s; 58°C 1 min; 72°C 2 min 30 s
			94°C 30 s; 58°C 1 min; 72°C 2 min 30 s
			94°C 30 s; 58°C 1 min; 72°C 2 min 30 s
第三轮	7	1	72°C 5 min
3 <sup>rd</sup> amplification			同第二轮(The same with secondary)