

表 2 SRAP 所用引物序列

Table 2 Primers sequence used in SRAP analysis

正向引物	序列	反向引物	序列
Forward primers	Sequence	Reverse primers	Sequence
me1	5'-TGAGTCCAAACCGGATA-3'	em1	5'-GACTGCGTACGAATTAAT-3'
me2	5'-TGAGTCCAAACCGGAGC-3'	em2	5'-GACTGCGTACGAATTTGC-3'
me3	5'-TGAGTCCAAACCGGAAT-3'	em3	5'-GACTGCGTACGAATTGAC-3'
me4	5'-TGAGTCCAAACCGGACC-3'	em4	5'-GACTGCGTACGAATTTGA-3'
me5	5'-TGAGTCCAAACCGGAAG-3'	em5	5'-GACTGCGTACGAATTAAC-3'
me6	5'-TGAGTCCAAACCGGTAA-3'	em6	5'-GACTGCGTACGAATTGCA-3'
me7	5'-TGAGTCCAAACCGGTCC-3'	em7	5'-GACTGCGTACGAATTCAA-3'
me8	5'-TGAGTCCAAACCGGTGC-3'	em8	5'-GACTGCGTACGAATTCTG-3'
me9	5'-TAGGTCCAAACCGGTAG-3'	em9	5'-GACTGCGTACGAATTCGA-3'
me10	5'-TAGGTCCAAACCGGTGT-3'	em10	5'-GACTGCGTACGAATTCAG-3'
		em11	5'-GACTGCGTACGAATTCCA-3'

