

表 5 RT-PCR 引物

Figure 5 Primers used for RT-PCR

	引物 Primer	引物序列(5'-3') Sequence (5'-3')	退火温度 (°C) AT (°C)	片段长度(bp) Product size (bp)	GenBank 登录号 GenBank accessaion No.
<i>Phlda2</i>	Phlda2-L	AGCTGGATAAGCGAAGTGAC	54	557	NM-001076521
	Phlda2-R	AGAAAAGAGAATGAACCGAG			
<i>GAPDH</i>	GAPDH-L	GCACAGTCAAGGCAGAGAAC	50	367	BTU85042
	GAPDH-R	GCGTGGACAGTGGTCATAAG			